ceramill® zirconia

 **ceramill® zi**

 **ceramill® zolid ht+**
white

 **ceramill® zolid ht+**
preshades

 **zolid gen-x**
multilayer

 **zolid drs**
multilayer

 **zolid bion**
multilayer

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EN Instruction Manual 28 - 51

FR Instructions d'utilisation 52 - 80

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AMANNGIRRBACH

- Original Gebrauchsanweisung -

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- Translation of the original Operating Instructions -

Table of Contents

Explanation of Symbols	29
General Safety Instructions	32
Suitable Personnel	32
Properties	33
Warranty/Exclusion of Liability	36
Testing the Blanks	36
Mounting to the Worktable	37
Application	37
Process steps in the dental laboratory	39
Process steps at the dentist's	48
Environmental Protection	51
Download information	51

EN

Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or to property damage are marked with the icon aside and are boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
–	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

Other symbols on the product

Symbol	Meaning
	Medical device
	Unique Device Identifier
	Catalog number
	Batch code
	Manufacturer
	Country of manufacture
	Consult instructions for use
	Use by YYYY-MM-DD or YYYY-MM
Rx only	Symbol that may be used in place of statement "CAUTION". US Federal law restricts this device to sale by or on the order of a dentist.
	Serial number
	Store in a dry place

General Safety Instructions

**CAUTION:**

Possible health impairment from zirconium-oxide dust!

- ▷ When processing, wear personal protective equipment (dust protection mask, safety glasses/goggles, ...).
- ▷ Do not breathe dust/fume/gas/mist/vapours/spray.
- ▷ Avoid release to the environment.

**NOTE:**

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Suitable Personnel

**NOTE:**

The product may only be processed by professionally trained dental technicians and dentists.

Properties

Intended Use

Zirconium-oxide blanks for permanent and partly removable dental prosthetics.

Product description

The blanks are made of zirconium oxide (ZrO_2) for type II, class 5 dental applications in accordance with DIN EN ISO 6872.

They are used for manufacturing permanent and partly removable prosthetic restorations (crowns and bridges, conical/ telescopic crowns, supra-constructions, and implant-supported restorations) using CNC milling machines (e.g. Ceramill milling machines).

After completion of the specified end-sintering, the final product meets the requirements of DIN EN ISO 6872.

Intended patient group: Suitable for patients of all ages and gender.

Clinical benefit of the product:

- _ Restoration of missing tooth and masticatory apparatus portions
- _ Restoration of chewing function and aesthetics

Safety data sheet/Declaration of conformity/SSCP

The safety data sheet and the SSCP of the products (according to the requirements of the Medical Device Regulation (MDR)) can be downloaded from the Amann Girrbach website under *Services > Downloads > Additional documents*. The declaration of conformity of the products are available upon request from the manufacturer.

Technical data

	Unit	Zi	Zolid Gen-X Multilayer
			Zolid DRS Multilayer
			Zolid HT+ Preshades
			Zolid Bion Multilayer
Flexural strength	MPa	≥ 800	≥ 800
E-module	GPa	≥ 200	≥ 200
Thermal expansion coefficient (CTE) (25 - 500 °C)	10 ⁻⁶ /K	10.4±0.5	10.5±0.5
Chemical solubility	µg/cm ²	< 100	< 100
Vickers hardness	HV10	1300±200	1300±200

The delivery scope can vary country-specific.

Chemical composition

Oxide	Zi	Mass percentage		
		Zolid Bion Multilayer	Zolid HT+ White	Zolid HT+ Preshadeds Zolid Gen-X Multilayer Zolid DRS Multilayer
ZrO ₂ + HfO ₂ + Y ₂ O ₃	≥ 99.0	≥ 99.0	≥ 99.0	≥ 99.0
Y ₂ O ₃	4.5 - 5.6	6.6 - 7.6	6.7 - 7.2	6.0 - 7.0
HfO ₂	≤ 5	≤ 5	≤ 5	≤ 5
Al ₂ O ₃	≤ 0.5	≤ 0.5	≤ 0.5	≤ 0.5
Other oxides	≤ 1	≤ 1	≤ 1	≤ 1

Shelf life of the blanks

When appropriately stored, the blanks can be used 5 years after the manufacturing date.

Storage

Store the blanks in the original packaging in a dry location.

Frequency of use

As an intermediate product, the respective product is not intended for single use. The final product (dental prosthesis) is patient-specific.

Warranty/Exclusion of Liability

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development and are tested according to the state of the art and to legal requirements. We reserve the right to make any resulting changes in handling and composition.

The respectively valid instruction-manual version can be found under www.amanngirrbach.com/instruction-manuals. This version replaces all previous versions.

Testing the Blanks

Upon receipt, it is imperative to visually check the proper condition of the blanks. After using a (transport-) damaged blank, complaint claims are no longer possible.

Mounting to the Worktable

The mounting of the blanks is described in the Operating Instructions of the respective milling machine.

Application

Indication ranges

- _ Anatomically reduced and fully anatomical (monolithic) crowns in the anterior and posterior tooth range (single-tooth crowns, inlays, onlays)
- _ Anatomically reduced and fully anatomical (monolithic) three to multi-unit bridges¹ with a maximum of three contiguous pontics in the anterior region and a maximum of two contiguous pontics in the posterior region
- _ Cantilever frames and bridges with a maximum of one bridge pontic (maximum one free-end pontic and no further than the second premolar)
- _ Implant-supported hybrid dentures (hybrid abutments and hybrid abutment crowns)
- _ Implant-supported, directly screwed restorations

Indication limitation for Canada:

- _ Single crowns
- _ Front-tooth bridges
- _ Posterior bridges with up to four units

1. Because of the Zolid DRS block form, only up to three-unit bridges can be fabricated with it.

Contraindications

- _ Insufficient tooth-structure availability
- _ Insufficient preparation
- _ Insufficient oral hygiene
- _ More than two connected bridge units in the posterior region, more than three connected intermediate units in the anterior region
- _ Known incompatibilities with respect to the components
- _ Heavily discoloured hard tooth structure
- _ Provisional insertion
- _ All indications that are not listed under “Indications”

Undesirable side effects/residual risks

- _ Mechanical failure (fracture of the restoration, chipping of the ceramic veneer or detachment of the restoration/loss of retention, loosening or breakage of the abutment screw and failure of the abutment or the suprastructure in the case of implant-supported care) with a low risk of potential swallowing or reversible damage to the mucous membranes.
- _ Biological complications (plaque accumulation, intolerance reactions, marginal/secondary caries, loosening of the supporting teeth) resulting in possible loss of the restoration.

Process steps in the dental laboratory

Material-specific frame parameters

The following material-specific frame parameters must be complied with when fabricating zirconium oxide frames in high-density sintered condition:

Minimum frame thickness in mm	Connector cross-section in mm ²		Maximum number of connected bridge units	
	Posterior area	Anterior area	Posterior area	Anterior area
0.5	≥ 7	≥ 9	3	2

Detailed minimum wall thicknesses and connector cross-sections

Indication	Maximum number of bridge units	Wall thickness in mm		Connector cross-section in mm ²
		Incisal/occlusal	Circular	
Primary components/Telescopic crown	–	0.7	0.5	–
Single crown	–	0.5	0.5	–
Front-tooth bridge	1	0.5	0.5	> 7
	3	0.7	0.5	> 9
Posterior bridge	1	0.7	0.5	> 9
	2	1.0	0.7	> 12
Cantilever bridge	1	1.0	0.7	> 12


Crown border parameters in the CAD software

▷ Pay attention to sufficient crown border thicknesses Recommended values:

- 1. Crown border thickness (horizontal): 0.2 mm
- 2. Angled: 0.3 mm
- 3. Angle: 60°
- 4. Height (vertical): 0 mm

For tangential preparations, increase the “crown border thickness” and “Angled” parameters, e.g.:


- 1. Crown border thickness (horizontal): 0.25 - 0.3 mm
- 2. Angled: 0.35 - 0.4 mm


Observe zirconium-oxide-efficient design forms and do not create pointed edges.

Positioning in the CAM software

- ▷ Apply enough sufficiently thick holding ridges (at least 3 ridges for a single crown; thickness >1.7 mm).


- ▷ Position holding ridges approx. 1 - 2 mm away from the crown margin, especially for long crown margins.

This provides support to the crown margin.

- ▷ Position the crown in the blank in such a manner that the clearance from the crown margin to the blank surface is at least 1 mm.

This ensures that the crown has no contact to the sinter skin of the blank.

Procedure when fabricating dental framework

For maximum fitting accuracy, the specific expansion/shrinkage value is given on the blanks.

When machining the blanks in milling machines, it is recommended to only use cutters with sufficient cutting edge quality.

- ▷ Depending on requirement in the respective CAM software, enter the appropriate value.
In case of uncertainties, contact your CAM supplier.
- ▷ Mill the blanks with the milling machine applying the milling strategies provided for the material.

Removal from the blank

Utmost care must be taken when removing the frames from the blank:

- ▷ Remove the frames only using suitable tools approved for zirconium oxide (avoid cutting discs).

For trimming and dressing the connectors and smoothing out the frames, use suitable tools approved for zirconium oxide.

For detailed information on the subject "Working zirconium oxide", please refer to the brochure "Processing Techniques Zolid". These can be found under www.amanngirrbach.com in the download area.

Colouration of the white zirconium oxide (Ceramill ZI, Ceramill Zolid HT+ White)

- ▷ Depending on requirement, partially or completely stain the white zirconium oxide frames prior to the sintering with dye solution (e.g. Ceramill Liquid new formula).

For detailed information on the subject “Colouration and finalising” please refer to the brochure “Processing Techniques Zolid”. These can be found under www.amanngirrbach.com in the download area.

Colouration of the pre-stained mono and polychrome zirconium oxide (Ceramill Zolid HT+ Preshade, Zolid Gen-X, Zolid DRS, Zolid Bion)

- ▷ Depending on requirement, partially stain the pre-colored mono and polychrome zirconia-oxide restorations (e.g. Ceramill Liquid new formula).
- ▷ After successful final sintering: color with stains or glaze materials or use as a framework for the veneering technique.

Final sintering

A high-temperature sintering furnace, e.g. Ceramill Therm or another high-quality and suitable furnace is recommended for the sintering (please also refer to the operating instructions of the respective sintering furnace).

For sintering Zolid DRS, use only the Ceramill Therm DRS high-temperature furnace.

Highly translucent zirconium oxide and opaque zirconium oxide can be sintered together in one sintering process.

NOTE:

When highly translucent zirconium oxide and opaque zirconium oxide are sintered together in the same sintering bowl and with the same sintering pearls, then this can lead to changes in the translucency. Therefore, it is advisable for highly translucent zirconium oxide and opaque zirconium oxide, to use individual sintering pearls in separate sintering bowls.

- ▷ Place the frames for sintering into the sintering bowl filled with sintering pearls.
- ▷ To ensure proper support, place the frames with slight pressure onto the sintering pearls.

NOTE:

- ▷ Pay attention that no pearls are jammed in the interdental spaces or located in the crown cavities!

▷ Perform the final sintering of the frames according to the following sintering parameters:

	Material	
	Ceramill Zi	
	Ceramill Zolid HT+ White	
	Ceramill Zolid HT+ Preshade	
Parameter	Zolid Gen-X	Zolid Bion
Heat-up phase	Room temperature to final temperature 1450 °C	Room temperature to final temperature 1500 °C
Heating rate	3 - 10 K/min	3 - 10 K/min
Dwell time at final temperature	2 hrs	2 hrs
Cooling phase	Final temperature 1450 °C to room temperature (at least < 200 °C); approx. 5 K/min (approx. 5 hours)	Final temperature 1500 °C to room temperature (at least < 200 °C); approx. 5 K/min (approx. 5 hours)

For detailed information on the subject “Sintering”, please refer to the brochure “Processing Techniques Zolid”. These can be found under www.amanngirrbach.com in the download area.

Frames for veneering

Post-processing

If required, the zirconium oxide frames can be reworked after the final sintering without pressure, using a water-cooled lab turbine and suitable diamond grinding points (recommended grain approx. 40 μm).

Frame preparation for veneering

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 50 μm , pressure ≤ 2 bar) and subsequent steaming off.

Veneer ceramics

For veneering with layer technique, a suitable veneering ceramic, which is matched to the coefficient of thermal expansion (CTE) of the framework material zirconium oxide, must be selected (see page 34). Suitable veneering ceramics for zirconium oxide include Creation ZI (Willy Geller) and Vita VM9 (Vita Zahnfabrik). Alternatively, zirconium oxide frameworks can be overpressed with veneering ceramic. Always observe the instructions for use of the zirconium dioxide material as well as the selected veneering material.

For bridges with more than five units and bridges with massive intermediate units, step-by-step long-term cooling to 500 °C is recommended.

Monolithic dentures

Post-processing

After the final sintering has taken place:

- ▷ If required, the zirconium oxide frames can be reworked without pressure, using a water-cooled lab turbine and suitable diamond grinding points (recommended grain approx. 40 µm).
- ▷ Polish the contact surfaces to the antagonist and to the neighbouring teeth to a high gloss when in high-density sintered condition.

Frame preparation for staining and glazing

After checking the fitting accuracy, contact points and occlusion:

- ▷ Clean the interior crown surfaces by means of corundum blasting (aluminium oxide, 50 µm, pressure ≤ 2 bar) and subsequent steaming off.

The exterior crown surfaces are not corundum blasted.

Staining and glazing

- ▷ Stain and glaze the units afterwards (e.g. with the Ceramill Stain & Glaze Kit).

For bridges with more than five units and bridges with massive intermediate units, step-by-step long-term cooling to 500 °C is recommended.

For detailed information on the subject “Staining and glazing” please refer to the brochure “Processing Techniques Zolid”. These can be found under www.amanngirrbach.com in the download area.

Process steps at the dentist's

Processing

All prosthetic components must be cleaned, disinfected and/or sterilised prior to each use (see “Guideline for Cleaning, Disinfection and Sterilization of Prosthetics Component Made of Amann Girrbach Materials”).

Fixation

Ceramill Zolid zirconium oxide restorations can be cemented conventionally with zinc phosphate cement or glass ionomer cement (e.g. GC Fuji Plus, GC) or adhesively using cement composites (Panavia[®] 21 or Panavia[®] F 2.0, Kuraray Noritake).

- ▷ With conventional cementing, ensure sufficient retention and a minimum stump height of 4 mm!

Temporary attachment is not recommended as scaffolding may be damaged upon removal.

Multilink Hybrid Abutment (Ivoclar Vivadent) is recommended for the extraoral attachment of a two-part abutment made of Ceramill Zolid zirconium oxide on a titanium adhesive base (hybrid abutment or a hybrid abutment crown). For attachment, follow the instructions of the manufacturer of the titanium adhesive base and the manufacturer of the attachment material.

Post-processing at the dentist's

When grinding is required by the dentist while installing the work into the patient's mouth, the units must be well re-polished to a high gloss. For grinding, we recommend using only diamond grinding points (recommended grain approx. 40 μm). For polishing, we recommend using diamond polishing bits.

NOTE:

When the work is not sufficiently polished, the antagonist may become damaged due to abrasion!

Follow-up examination at the dentist's

It is urgently recommended to check the dentures in the patient's mouth once a year. In this, the remaining dentition, the antagonists and the soft tissue are also to be inspected. Corrective measures are to be undertaken as required. Here, it is also important that the objects are polished afterwards to a high gloss.

Information to be passed on to the patient

(by healthcare professionals).

Possible complications and residual risks for zirconium-oxide-based restorations in general and not specific to Amann Girrbach restorative materials are:

- _ Mechanical failure (fracture of the restoration, chipping of the ceramic veneer or detachment of the restoration/loss of retention, loosening or breakage of the abutment screw and failure of the abutment or the suprastructure in the case of implant-supported care) with a low risk of potential swallowing or reversible damage to the mucous membranes.
- _ Biological complications (plaque accumulation, intolerance reactions, marginal/secondary caries, loosening of the supporting teeth) resulting in possible loss of the restoration.

To reduce possible complications/residual risks, it is recommended to have the denture checked by the dentist during the regular dental check-up (once a year).

Environmental Protection

Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally-friendly and recyclable.

Disposal of the blanks

Larger quantities may not be disposed of together with household waste. Smaller quantities may be disposed of together with household waste. Do not dispose of into sewer systems. Dispose of non-cleanable packaging properly.

Always dispose of according to official regulations.

Download information

Further instructions/manuals can be downloaded under www.amanngirrbach.com/instruction-manuals.

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Amann Girrbach AG
Amann Girrbach GmbH