DE	Gebrauchsanweisung	3
EN	Instruction Manual	9
FR	Instructions d'utilisation	15
IT	Istruzioni per l'uso	21
ES	Instrucciones de uso	27
CS	Návod k použití	33
EL	Οδηγίες χειρισμού	38
HR	Upute za uporabu	44
HU	Használati utasítás	50
PL	Instrukcja obsługi	56
PT	Instruções de uso	62
RO	Instrucţiuni de utilizare	68
SV	Bruksanvisning	74
DA	Brugsanvisning	80
ET	Kasutusjuhend	86
LT	Naudojimo instrukcija	92
LV	Lietošanas pamācība	98
NL	Gebruiksaanwijzing	104
SK	Návod na používanie	110
SL	Navodila za uporabo	116
TR	Kullanıcı kılavuzu	122

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Instruction Manual

PRODUCT NAME	Ceramill® CoCr
DESCRIPTION	NPM Cobalt-based dental milling alloy for metal-ceramic, type 4
DIMENSION	Ø 98.3 mm: 10 mm, 12 mm, 13.5 mm, 16 mm, 18 mm, 20 mm
CONTENT	1 Disc

CHEMICAL COMPOSITION (typical values):

Co %	Cr %	W %	Si %	Mn/Fe %
61.7	27.8	8.5	1.6	≤ 1.0

TYPICAL TECHNICAL DATA:

Yield strength 0.2 %	375 MPa
Elongation	14,0 %
Tensile strength	550 MPa
E-module	210 GPa
Density	8.6 g/cm ³
Corrosion resistance	≤ 200 µg/cm ² /7d
Hardness	285 HV 10/30
CTE (25 - 500 °C)	14.5 x 10 ⁻⁶ K ⁻¹
Max. firing temp.	980 °C
APPLIED NORM:	DIN EN ISO 22674
	ED GmbH is certified according to DIN EN ISO 13485

Intended use

Ceramill® CoCr is a medical device for the milling of crowns and bridges.

Only for professional user (Dental Technician, Dentist)!

The intended patient group provides for persons with partially or non-dentate jaws.

Indication

Crowns and bridges for metal-ceramics up to 16 units (max. 4 units in the anterior region and max. 3 units in the posterior region), superstructures, bars.

Contraindication

In case of known allergic reactions to any of the ingredients.

Milling

Ceramill® CoCr is construed for CNC milling machines. Please follow the instructions and parameters of the respective manufacturer of CAM Software and the CNC milling machine.

CAD

The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Choose a sufficient connector dimension (6 - 9 mm²). Sharp edges and undercuts should be avoided.

Cutting out the frameworks from the Blank

Remove the milled frameworks with suitable cutting tools and smoothing the supports.

Soldering/Laser welding

Ceramill® CoCr can be soldered with all suitable solder.

Ceramill® CoCr parts should not be soldered with gold or palladium solder. **Ceramill® CoCr** is also ideally suitable for laser-welding.

Preparation before ceramic veneering

The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with minimum 110 µm of Aluminium oxide with 3-4 bar and clean with steam cleaner. Oxide firing is not mandatory but can be done as an option for 5 minutes at 980 °C (cleaning firing). The frame needs to be sandblasted again to remove the present oxide layer thoroughly. In the end the cleaning by steam cleaner is mandatory. If you use a ceramic bonder please consider the instruction for use of the manufacturer.

Handling conditions/Safety

Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 - DIN EN 149.

Residual risks and side effects

If the instructions are observed during the production processes, incompatibilities with CoCr alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

Disinfection of the dental prosthesis before insertion

Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Single-use

Used milling disc should be not processed for further use as a medical device.

Disposal Instructions

Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions

Temperature, humidity or light has no effect on the product properties

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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Amann Girrbach AG