# CLEARFIL™ Universal Bond Quick

<b>EN</b>	ENGLISH	English	DENTAL BONDING AGENT
<b>FR</b>	FRANÇAIS	French	AGENT DE LIAISON DENTAIRE
<b>ES</b>	ESPAÑOL	Spanish	AGENTE DE ADHESIÓN DE USO DENTAL
<b>IT</b>	ITALIANO	Italian	AGENTE DI ADESIONE DENTALE
<b>NL</b>	NEDERLANDS	Dutch	DENTALE BONDING
<b>DE</b>	DEUTSCH	German	Dentaler Haftvermittler
<b>SV</b>	SVENSKA	Swedish	DENTALT BONDINGSMATERIAL
<b>NO</b>	NORSK	Norwegian	DENTALT BINDINGSMIDDEL
<b>FI</b>	SUOMI	Finnish	HAMPAAN SIDOSAINE
<b>DA</b>	DANSK	Danish	DENTALT BONDINGMATERIALE
<b>PT</b>	PORTUGUÊS	Portuguese	AGENTE DE ADESÃO DENTÁRIO
<b>EL</b>	ΕΛΛΗΝΙΚΑ	Greek	ΟΔΟΝΤΙΚΟΣ ΣΥΓΚΟΛΛΗΤΙΚΟΣ ΠΑΡΑΓΟΝΤΑΣ
<b>TR</b>	TÜRKÇE	Turkish	DENTAL BONDING AJAN
<b>PL</b>	POLSKI	Polish	STOMATOLOGICZNY SYSTEM WIĄŻĄCY
<b>RO</b>	ROMÂNĂ	Romanian	AGENT DENTAR ADEZIV
<b>HR</b>	HRVATSKI	Croatian	DENTALNI ADHEZIV
<b>HU</b>	MAGYAR	Hungarian	FOGÁSZATI KÖTŐANYAG
<b>SL</b>	SLOVENŠČINA	Slovenian	DENTALNO VEZIVNO SREDSTVO
<b>CS</b>	ČESKY	Czech	DENTÁLNÍ POJIVO
<b>SK</b>	SLOVENSKY	Slovak	DENTÁLNE TMELIACE ČINIDLO
<b>BG</b>	БЪЛГАРСКИ	Bulgarian	СТОМАТОЛОГИЧЕН БОНДИНГ АГЕНТ
<b>UK</b>	УКРАЇНСЬКА	Ukrainian	СТОМАТОЛОГІЧНИЙ АДГЕЗИВ
<b>ET</b>	EESTI	Estonian	DENTAALNE ADHESIIV
<b>LV</b>	LATVIEŠU	Latvian	ZOBĀRSTNIECĪBAS SAISTVIELA
<b>LT</b>	LIETUVIŠKAI	Lithuanian	ODONTOLOGINIS RIŠIKLIS

<b>MD</b>	<b>EN</b>	Medical Device	<b>DA</b>	Medicinsk udstyr	<b>CS</b>	Zdravotnický prostředek
	<b>FR</b>	Dispositif médical	<b>PT</b>	Dispositivo médico	<b>SK</b>	Zdravotnícka pomôcka
	<b>ES</b>	Producto sanitario	<b>EL</b>	Ιατροτεχνολογικό προϊόν	<b>BG</b>	Медицинско изделие
	<b>IT</b>	Dispositivo medico	<b>TR</b>	Tıbbi cihaz	<b>UK</b>	Виріб медичного призначення
	<b>NL</b>	Medisch hulpmiddel	<b>PL</b>	Wyrób medyczny	<b>ET</b>	Meditsiiniseade
	<b>DE</b>	Medizinprodukt	<b>RO</b>	Dispozitiv medical	<b>LV</b>	Medicīniska ierīce
	<b>SV</b>	Medicinteknisk produkt	<b>HR</b>	Medicinski proizvod	<b>LT</b>	Medicinos priemonė
	<b>NO</b>	Medisinsk utstyr	<b>HU</b>	Orvostechikai eszköz		
	<b>FI</b>	Lääkinnällinen laite	<b>SL</b>	Medicinski pripomoček		

**Kuraray Noritake Dental Inc.**

1621 Sakazu, Kurashiki, Okayama 710-0801, Japan

**Kuraray Europe GmbH** (Importer)

Philipp-Reis-Str. 4,


65795 Hattersheim am Main, Germany

Phone:+49 (0)69 305 35835 Fax:+49 (0)69 305 98 35835

URL:<https://www.kuraraynoritake.eu>

1562R768R

**CLEARFIL™ Universal Bond Quick****I. INTRODUCTION**

"CLEARFIL Universal Bond Quick" consists of BOND and K-ETCHANT Syringe. The BOND is a light curing bonding agent that allows the treatment of dentin, enamel, and prosthetic materials. Depending on the indication, the BOND is used as self-etching or with K-ETCHANT Syringe for selective enamel etching or total-etching procedures. The BOND is intended to be used for both direct and indirect restorations. "CLEARFIL DC Activator" activates the dual-curing mechanism of the BOND; however, the addition of "CLEARFIL DC Activator" to the adhesive is not required when using it with "CLEARFIL DC CORE PLUS" or "PANAVIA SA Cement Plus". The BOND is available in both Bottle and Unit Dose delivery systems. K-ETCHANT Syringe is an etching gel that consists of 35 % phosphoric acid aqueous solution and colloidal silica.

**II. INDICATIONS**

"CLEARFIL Universal Bond Quick" is indicated for the following uses:

- [1] Direct restorations using light-cured composite resin
- [2] Cavity sealing as a pretreatment for indirect restorations
- [3] Treatment of exposed root surfaces
- [4] Treatment of hypersensitive teeth
- [5] Intraoral repairs of fractured restorations
- [6] Post cementation and core build-ups
- [7] Cementation of indirect restorations

**III. CONTRAINDICATIONS**

Patients with a history of hypersensitivity to methacrylate monomers and this product

**IV. POSSIBLE SIDE EFFECTS**

- [1] The oral mucosal membrane may turn whitish when contacted by the BOND due to the coagulation of protein. This is usually a temporary phenomenon that will disappear in a few days. Instruct patients to avoid irritating the affected area while brushing.
- [2] K-ETCHANT Syringe may cause inflammation or erosion due to its chemistry.

**V. INCOMPATIBILITIES**

- [1] Do not use eugenol-containing materials for pulp protection or temporary sealing, since the eugenol can retard the curing process.
- [2] Do not use hemostatics containing ferric compounds, since these materials may impair adhesion and may cause discoloration of the tooth margin or surrounding gingiva, due to remaining ferric ions.
- [3] When using hemostatics containing aluminum chloride, minimize the quantity; use caution to prevent contact with the adherent surface. Failure to do so might weaken the bond strength to the tooth structure.

**VI. PRECAUTIONS****1. Safety precautions**

1. This product contains substances that may cause allergic reactions. Avoid use of the product in patients with known allergies to methacrylate monomers or any other components in the product.
2. If the patient demonstrates a hypersensitivity reaction such as rash, eczema, features of inflammation, ulcer, swelling, itching, or numbness, discontinue use of the product, remove the product and seek medical attention.
3. Use caution to prevent the product from coming in contact with the skin or getting into patient's eyes. Before using the product, cover the patient's eyes with a towel to protect them in the event of splashing material.
4. If the product comes in contact with human body tissues, take the following actions:
  - <If the product gets in the eye>  
Immediately wash the eye with copious amounts of water and consult a physician.
  - <If the product comes in contact with the skin or the oral mucosa>  
Immediately wipe the area with a cotton pellet or a gauze pad moistened with alcohol and rinse with copious amounts of water.
5. Use caution to prevent the patient from accidentally swallowing the product.
6. Avoid looking directly at the dental curing light when curing the BOND.
7. Avoid using the same BOND dispensed into a well of the dispensing dish, Unit Dose and applicator brush for different patients to prevent cross contamination. Unit Dose and applicator brushes are single use only. Discard them after use. The needle tip is single-use only. Do not reuse it to prevent cross-contamination. Discard it after use.
8. Wear gloves or take other appropriate protective measures to prevent the occurrence of hypersensitivity that may result from contact with methacrylate monomers or any other components in the product.
9. If the instruments associated with this product are damaged, use caution and protect yourself; immediately discontinue use.
10. Dispose of this product as a medical waste to prevent infection. The needle tip must be disposed of after covering the tip of the needle to prevent injury.
11. If a dental temporary filling/temporary cement is a resin-based material which is applied on the surface of cavity sealing with BOND or coating with a composite resin, use a dental separation material according to the manufacturer's instructions to avoid bonding between the temporary material and the surface.

**2. Handling and manipulation precautions****[Common precautions]**

1. This product must not be used for purposes other than specified in [II.INDICATIONS].
2. The use of this product is restricted to licensed dental professionals.
3. Use a pulp capping agent in a cavity close to the pulp or in the event of accidental pulp exposure.

**[BOND]**

1. The BOND contains ethanol, a flammable substance. Do not use it near an open flame.
2. To prevent poor performance or poor handling characteristics, observe the specified light-curing times and other handling requirements.
3. Clean the cavity sufficiently to prevent poor bonding. If the adherent surface is contaminated with saliva or blood, wash it thoroughly and dry before bonding.
4. **Bottle:** Use the light-blocking plate to avoid exposing the material to an operating light or ambient light, and use within the working time indicated below after dispensing or mixing.

Material	Working time
BOND	7 minutes
BOND + CLEARFIL DC Activator	90 seconds

As the ethanol in BOND evaporates, the viscosity increases, thereby making it difficult to apply.

**Unit Dose:** Apply BOND with the applicator brush immediately after snapping off the container cap.

5. During the application of BOND to the adherent surface with rubbing motion, move the lighting-spot out of the mouth or turn off the light to prevent the applied BOND from being exposed to the operating light. In addition, after application of BOND, dry sufficiently until BOND does not move by blowing mild air.
6. BOND contains ethanol and water. "CLEARFIL DC Activator" contains ethanol. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until BOND or the mixture of BOND and "CLEARFIL DC Activator" does not move; otherwise, the adhesion effect will be impaired. In order to dry sufficiently, adjust the air pressure according to the shape and size of the cavity and the prosthetic appliance. Use a vacuum aspirator to prevent BOND or the mixture from scattering.
7. If the treated surface is contaminated, wash it with water, dry, or clean with alcohol, and treat again with BOND.
8. **Bottle:** Do not mix BOND with other bonding agents except "CLEARFIL DC Activator".  
**Unit Dose:** Do not mix BOND with other bonding agents. Do not use Unit Dose with CLEARFIL DC Activator due to the design of the Unit Dose container.
9. **Unit Dose:** When using for core build-up or cementation, use BOND only with "CLEARFIL DC CORE PLUS" or "PANAVIA SA CEMENT Plus".
10. **Bottle:** Light-cure the mixture of the BOND and "CLEARFIL DC Activator", otherwise, the working time will be shortened dramatically.
11. **Bottle:** The container should be tightly capped immediately after use to reduce the evaporation of the volatile solvent (ethanol contained in the BOND). If the liquid does not flow easily from the nozzle, do not dispense the BOND forcibly from the plugged container.
12. **Bottle:** If the BOND has not been used in a long time, the BOND may not flow easily; shake the container before use.

**[K-ETCHANT Syringe]**

1. Be careful not to contaminate it with saliva or blood. If the treated surface is contaminated, re-treat.
2. Be careful to avoid cross-contamination. Disinfect the syringe by wiping it with an absorbent cotton with alcohol both before and after use. Cover the entire syringe with a disposable plastic barrier to prevent saliva and blood contamination.
3. If the product adheres to clothing, wash it off with water.
4. After each use, remove the needle tip from the syringe and recap the syringe immediately and tightly.
5. Etching the vital dentin may cause post-operative sensitivity.

**[Dental light-curing unit]**

1. Low light intensity causes poor adhesion. Check the lamp for service life and the dental curing light guide tip for contamination. It is advisable to check the dental curing light intensity using an appropriate light evaluating device at appropriate intervals.
2. The emitting tip of the dental curing unit should be held as near and vertical to the resin surface as possible. If a large resin surface has to be light-cured, it is advisable to divide the area into several sections and light-cure each section separately.

**3. Storage precautions**

1. The product must be used by the expiration date indicated on the package.
2. BOND must be refrigerated (2-8°C/ 36-46°F) when not in use, and should be brought to room temperature for more than 15 minutes before using. Especially, Bottle must be left standing until it comes to room temperature after taking out of the refrigerator; otherwise, an excessive amount of liquid might be dispensed or the liquid may ooze after use.
3. K-ETCHANT Syringe must be stored at 2 - 25°C / 36 - 77°F when not in use.
4. Keep away from extreme heat, direct sunlight or a flame.
5. The product must be stored in a proper place where only dental practitioners can access.

**VII. COMPONENTS**

Please see the outside of the package for contents and quantity.

<Principal ingredients>

- 1) BOND
  - 10-Methacryloyloxydecyl dihydrogen phosphate (MDP)
  - Bisphenol A diglycidylmethacrylate (Bis-GMA)
  - 2-Hydroxyethyl methacrylate (HEMA)
  - Hydrophilic amide monomers
  - Colloidal silica
  - Silane coupling agent
  - Sodium fluoride
  - dl-Camphorquinone
  - Ethanol
  - Water
- 2) K-ETCHANT Syringe
  - Phosphoric acid
  - Water
  - Colloidal Silica
  - Pigment

- 3) Accessories
- Applicator brush (fine <silver>)
  - Dispensing dish
  - Light-blocking plate
  - Needle tip (E)

## VIII. CLINICAL PROCEDURES

### A. Standard procedure I

- [1] Direct restorations using light-cured composite resin  
 [2] Cavity sealing as a pretreatment for indirect restorations  
 [3] Treatment of exposed root surfaces  
 [4] Treatment of hypersensitive teeth

#### A-1. Isolation and Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

#### A-2. Cavity or root surface preparations

Remove any infected dentin and prepare the cavity in the usual manner. When treating hypersensitive teeth, clean the root surface as usual. Clean thoroughly by spraying with water, and dry with air or with cotton pellets.

#### A-3. Pulp protection

Any actual or near pulp exposure can be covered with a hard setting calcium hydroxide material. There is no need for cement lining or base. Do not use eugenol materials for pulp protection.

#### A-4. Pretreatment of tooth

Choose one of the three etching procedures before applying BOND.

[NOTE]

When treating hypersensitive teeth, select section A-4a before applying BOND.

##### A-4a. Self-etching procedure

Move to section A-5 without etching with K-ETCHANT Syringe.

##### A-4b. Selective enamel etching procedure

Apply K-ETCHANT Syringe to the uncut and/or cut enamel. Leave it in place for 10 seconds, then rinse and dry.

##### A-4c. Total-etching procedure

Apply K-ETCHANT Syringe to the entire cavity (enamel and dentin), leave it in place for 10 seconds, then rinse and dry.

#### A-5. Application of BOND

1. **Bottle:** Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.

[CAUTION]

Use the light-blocking plate to avoid exposing the material to an operating light or ambient light; use within 7 minutes after dispensing.

**Unit Dose:** Snap off the container cap.

[CAUTION]

When snapping off the container cap, do not tilt in order to avoid spilling BOND.

2. Apply BOND with a rubbing motion to the entire cavity wall with the applicator brush. No waiting time is required.

[NOTE]

Use caution not to allow saliva or exudate to contact the treated surfaces.

3. Dry the entire cavity wall sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

4. Light-cure BOND with a dental curing unit. (See table "Dental curing unit and curing time").

Table: Dental curing unit and curing time

Type	Light source	Light Intensity	Light-curing time
Halogen	Halogen lamp	More than 400 mW/cm <sup>2</sup>	10 seconds
LED	BLUE LED*	800 – 1400 mW/cm <sup>2</sup>	10 seconds
		More than 1500 mW/cm <sup>2</sup>	5 seconds

The effective wavelength range of each dental curing unit must be 400 - 515 nm.

\*Peak of emission spectrum: 450 - 480 nm

#### A-6. Placement of composite resin restorative material; treatment of hypersensitive teeth; or cavity sealing.

##### A-6a. Direct restorations using light-cured composite resin.

Place composite resin (e.g. CLEARFIL MAJESTY ES-2, CLEARFIL MAJESTY ES Flow) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

##### A-6b. Treatment of exposed root surfaces

Place a thin coat of composite resin (e.g. CLEARFIL MAJESTY ES Flow) onto the tooth, light-cure according to the manufacturer's instructions. Remove un-polymerized resin with a cotton pellet or a gauze moistened with alcohol.

##### A-6c. Treatment of hypersensitive teeth

Remove the un-polymerized layer of BOND with a cotton pellet or a gauze moistened with alcohol.

##### A-6d. Cavity sealing

If necessary, place a thin coat of composite resin (e.g. CLEARFIL MAJESTY ES Flow) onto the tooth, light-cure according to the manufacturer's instructions.

Remove un-polymerized resin of the bonding or the composite resin surface with a cotton pellet or a gauze moistened with alcohol.

### B. Standard procedure II

#### [5] Intraoral repair of fractured restorations

##### B-1. Preparation of adherent surfaces

Roughen the adherent surfaces using a diamond point or blasting with 30 to 50 µm alumina powder at air pressure of 0.1-0.4 MPa (14-58 PSI/ 1-4 kgf/cm<sup>2</sup>). The air pressure should be properly adjusted to suit the material and/or shape of the prosthetic restoration, using caution to prevent it from chipping. Place a bevel at the marginal area.

##### B-2. K-ETCHANT Syringe treating of adherent surfaces

Apply K-ETCHANT Syringe to the adherent surface (including tooth structure). Leave in place for 5 seconds, then rinse and dry.

##### B-3. Application of BOND

Apply BOND with a rubbing motion to the entire adherent surface. Refer to section A-5.

[NOTE]

For optimal performance, apply a silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER PLUS) to the surface of the silica-based glass ceramic (e.g. conventional porcelain, Lithium disilicate) and apply a metal adhesive primer (e.g. ALLOY PRIMER) to the surface including precious metal according to the manufacturer's instructions BEFORE applying BOND.

##### B-4. Place composite resin restorative

Place composite resin (e.g. CLEARFIL MAJESTY ES-2) into the cavity, light-cure, finish and polish according to the manufacturer's instructions.

[NOTE]

Use an opaque resin (e.g. CLEARFIL ST OPAQUER) before composite resin placement to mask metal color.

### C. Standard procedure III

#### [6] Post cementation and core build-ups

When using with "CLEARFIL DC CORE PLUS", the use of "CLEARFIL DC Activator" is not necessary.

##### C-1. Isolation and Moisture control

Avoid contamination of the treatment area from saliva or blood to produce optimal results. A rubber dam is recommended to keep the tooth clean and dry.

##### C-2. Preparing the root canal

Prepare and clean the root canal in the usual manner.

##### C-3. Post preparation

Choose either C-3a or C-3b based on the post used. Please follow the Instructions for Use of the restorative material. In the absence of specific instructions, we recommend the following procedure:

##### C-3a. For Glass Fiber Posts

Apply K-ETCHANT Syringe to the post surface. Leave it in place for 5 seconds, then rinse and dry.

[CAUTION]

- Do not blast glass fiber posts with alumina powder, or the posts could be damaged.
- Avoid all contamination of the surfaces to be treated during pre-treatment and until the final core build-up.

##### C-3b. For Metal Posts

Roughen the adherent surface by blasting with 30 to 50 µm alumina powder at an air pressure of 0.2-0.4 MPa (29-58 PSI/ 2-4 kgf/cm<sup>2</sup>). The air pressure should be properly adjusted to suit the material. After blasting, clean the prosthetic restoration by using ultrasonic cleaner for 2 minutes followed by drying it with an air stream.

##### C-4. Post surface treatment

Choose either procedure based on the material used.

##### C-4a. When using with CLEARFIL DC CORE PLUS

1. Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.
2. Apply BOND to the entire post surface with the applicator brush.
3. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

##### C-4b. When using with other dual-/self-cured resin core materials or dual-/self-cured resin cements

1. **Bottle:** Dispense one drop each of BOND and "CLEARFIL DC Activator" into a well of the dispensing dish and mix them with the applicator brush.

[CAUTION]

Use the light-blocking plate to avoid exposing the material to an operating light or ambient light; use within 90 seconds after mixing.

**Unit Dose:** Not applicable for this use.

2. Apply the mixture to the post surface.

3. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent the mixture from scattering.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

4. Light-cure the mixture with a dental curing unit (See table "Dental curing unit and curing time").

[CAUTION]

Working time will be shortened dramatically when not light-curing the mixture on the post.

##### C-5. Pretreatment of tooth

Apply K-ETCHANT Syringe, as needed. Refer to section A-4.

##### C-6. Bonding

Choose either procedure based on the material used.

##### C-6a. When using with "CLEARFIL DC CORE PLUS"

1. Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.
2. Apply BOND with a rubbing motion to the entire cavity wall with the applicator brush. No waiting time is required.

[NOTE]

Use caution not to allow saliva or exudate to contact the treated surfaces.

3. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering. Remove the excess BOND with a paper point. After wiping away the excess BOND, dry the adherent surface again if necessary.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

4. Light-cure BOND with a dental curing unit (See table "Dental curing unit and curing time").

#### **C-6b. When using with other dual-/self-cured resin core materials or dual-/self-cured resin cements**

1. Dispense one drop each of BOND and "CLEARFIL DC Activator" into a well of the dispensing dish and mix them with the applicator brush.

[CAUTION]

Use the light-blocking plate to avoid exposing the material to an operating light or ambient light; use within 90 seconds after mixing.

2. Apply the mixture with a rubbing motion to the entire cavity wall with the applicator brush. No waiting time is required.

[NOTE]

Use caution not to allow saliva or exudate to contact the treated surfaces.

3. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent the mixture from scattering. Remove the excess mixture with a paper point. After wiping away the excess mixture, dry the adherent surface again if necessary.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

4. Light-cure the mixture with a dental curing unit (See table "Dental curing unit and curing time").

[CAUTION]

Working time will be shortened dramatically when not light-curing the mixture in the root canal.

#### **C-7. Placing the post and core build-up**

Place the post and core build-up using "CLEARFIL DC CORE PLUS", or other resin material according to the manufacturer's instructions.

#### **D. Standard procedure IV**

##### **[7] Cementation of indirect restorations**

When using with "PANAVIA SA Cement Plus", the use of "CLEARFIL DC Activator" is not necessary.

##### **D-1. Conditioning the cavity and core (tooth, metal, composite) surfaces**

1. Remove the temporary sealing material and temporary cement in the usual manner, and clean the cavity with moisture control.
2. Trial fit the prosthetic restoration to check its fit on the cavity or core (tooth, metal, composite). When using a try-in paste to check the color, follow the manufacturer's instructions.

##### **D-2. Surface preparation of prosthetic restorations**

Choose either D-2a or D-2b based on the restoration used. Follow the Instructions for Use of the restorative material. In the absence of specific instructions, we recommend the following procedure:

###### **D-2a. For silica-based glass ceramic (e.g. conventional porcelain, Lithium disilicate)**

Etch the glass ceramic surfaces with hydrofluoric acid in accordance with the manufacturer's instructions, and thoroughly wash and dry the surface.

###### **D-2b. For Metal oxides (e.g. Zirconia), metals or composite resins**

Roughen the adherent surface by blasting with 30 to 50µm alumina powder at an air pressure of 0.2-0.4 MPa (29-58 PSI/ 2-4 kgf/cm<sup>2</sup>). The air pressure should be properly adjusted to suit the material and/or shape of the prosthetic restoration, using caution to prevent chipping. After blasting, clean the prosthetic restoration by using an ultrasonic cleaner for 2 minutes, followed by drying it with an air stream.

##### **D-3. Pretreatment of prosthetic restorations**

Choose either procedure based on the material used.

###### **D-3a. When using with PANAVIA SA Cement Plus**

1. Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.
2. Apply BOND to the entire adherent surface with the applicator brush.
3. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE]

- To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

- For optimal performance, apply a silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER PLUS) to the surface of the silica-based glass ceramic (e.g. conventional porcelain, Lithium disilicate) instead of BOND, according to the manufacturer's instructions.

- Application of BOND to metal oxides or metals is not necessary since "PANAVIA SA Cement Plus" bonds strongly to these surfaces.

###### **D-3b. When using with other dual- or self-cured resin cements**

1. **Bottle:** Dispense one drop each of BOND and "CLEARFIL DC Activator" into a well of the dispensing dish and mix them with the applicator brush.

[CAUTION]

Use the light-blocking plate to avoid exposing the material to an operating light or ambient light; use within 90 seconds after mixing.

**Unit Dose:** Not applicable for this use.

2. Apply the mixture to the adherent surface.
3. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent the mixture from scattering.

[NOTE]

- To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

- For optimal performance, apply a silane coupling agent (e.g. CLEARFIL CERAMIC PRIMER PLUS) to the surface of the silica-based glass ceramic (e.g. conventional porcelain, Lithium disilicate) instead of the mixture, according to the manufacturer's instructions.

4. Light-cure the mixture with a dental curing unit (See table "Dental curing unit and curing time").

[CAUTION]

Working time will be shortened dramatically when not light-curing the mixture on the adherent surface.

#### **D-4. Pretreatment of tooth**

As necessary, apply K-ETCHANT Syringe. Refer to section A-4.

#### **D-5. Bonding**

Choose either procedure based on the material used.

##### **D-5a. When using with PANAVIA SA Cement Plus**

1. Dispense the necessary amount of BOND into a well of the dispensing dish immediately before application.
2. Apply BOND with a rubbing motion to the entire cavity wall with the applicator brush. No waiting time is required.

[NOTE]

Use caution not to allow saliva or exudate to contact the treated surfaces.

3. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until BOND does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

##### **D-5b. When using with other dual- or self-cured resin cement**

1. Dispense one drop each of BOND and "CLEARFIL DC Activator" into a well of the dispensing dish and mix them with the applicator brush.

[CAUTION]

Use the light-blocking plate to avoid exposing the material to an operating light or ambient light; use within 90 seconds after mixing.

2. Apply the mixture with a rubbing motion to the entire cavity wall with the applicator brush. No waiting time is required.

[NOTE]

Use caution not to allow saliva or exudate to contact the treated surfaces.

3. Dry the entire adherent surface sufficiently by blowing mild air for more than 5 seconds until the mixture does not move. Use a vacuum aspirator to prevent BOND from scattering.

[NOTE]

To dry thoroughly, adjust the air pressure according to the shape and size of the adherent surface.

4. Light-cure the mixture with a dental curing unit (See table "Dental curing unit and curing time").

[CAUTION]

Working time will be shortened dramatically when not light-curing the mixture on the adherent surface.

#### **D-6. Cementing**

Cement the prosthetic restoration using "PANAVIA SA Cement Plus" or other resin cement according to the manufacturer's instructions.

[NOTE]

When using a partial light-curing (or "Tack-Cure") technique, the setting time of the excess cement will be faster.

The BOND or the mixture of BOND with CLEARFIL DC Activator may accelerate the light-curing of the cement.

[WARRANTY]

Kuraray Noritake Dental Inc. will replace any product that is proven to be defective. Kuraray Noritake Dental Inc. does not accept liability for any loss or damage, direct, consequential or special, arising out of the application or use of or the inability to use these products. Before using, the user shall determine the suitability of the products for the intended use and the user assumes all risk and liability whatsoever in connection therewith.

[NOTE]


If a serious accident attributable to this product occurs, report it to the manufacturer's authorized representative shown below and the regulatory authorities of the country in which the user/patient resides.

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"CLEARFIL", "CLEARFIL MAJESTY", "CLEARFIL ST", "CLEARFIL DC CORE PLUS" and "PANAVIA" are trademarks of KURARAY CO., LTD.