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Crown-Click de Luxe is a mechanical crown and bridge remover for the removal of fixed and temporary dentures.

- Using an intensive impulse of the impact stress inside the instrument, the cement layer breaks and the denture can be removed easily. The periodontium does not get damaged in this process.
- Easy handling: One hand holds the instrument, the other hand fixates the hook top at the denture with practically no assistance needed.
- The infinitely adjustable impact intensity as well as three different hooks guarantee a universal use.

Function: The impact strength is adjusted by turning the knurled-head screw (visible in the empty space of the instrument shelter).

- Stronger: Turn the knurled-head screw towards the hook – raises spring tension
- Weaker: Turn the knurled-head screw towards the trigger – releases spring tension

Changing instruments: The hooks can be exchanged easily and fast. Put the delivered fork wrench on the groove of the hook fixing aid and penetrate into the drill of the receiver hook with the guide pin. By holding the fixing tube with the fork wrench and turning the guide pin to the left, the concerned hook can be taken out easily. The complete device is made of high-quality stainless steel and practically maintenance-free.

Scope of supply:

Set: Device with all three instruments (incl. fork wrench, guide pin for mounting and demounting the hooks)

Tip no. 1 long

Tip no. 2 short

Tip no. 3 short and angles

REF 452 205

REF 452 206

REF 452 207

REF 452 208

Processing instructions according to EN ISO 17664:201

Warnings: Attention! The product must be cleaned, disinfected and sterilized according to these processing in-

structions before its first use as well as after every subsequent use. Disinfection alone is not sufficient. Thorough cleaning and disinfection are decisive prerequisites for effective sterilization. The processing should be started as soon as possible but not later than 2 hours after use. During the processing, the product should not be unnecessarily exposed to moisture or humidity. Please also observe the legal regulations applicable in your country as well as the hygiene regulations of the medical practice and/or the hospital. **Limitations in the processing:** Frequent processing has only minor effects on this product. The end of the product's useful life is mainly determined by wear and damage caused by use. In case of doubt, the products should always be sorted out and replaced early. The user alone is responsible for making the decision regarding repeated use of the product. If the product is used too frequently, the manufacturer does not accept any warranty for the function, performance and safety of the product.

Personal protective equipment: For reasons of occupational safety and to minimize infections/cross infections, suitable personal protective equipment (protective clothes, protective gloves, protective goggles and nose/mouth protection) must be used during the entire processing.

INSTRUCTIONS:

INITIAL TREATMENT AT THE PLACE OF USE

Remove surface contamination using a disposable cloth/paper cloth when you are still at the place of use. Rinse the products with water (at least drinking water quality) no later than 2 hours after use. Do not allow residues or contamination of any kind to dry on the product. When using the product for the first time at the place of use, do not use any aldehyde-containing or alcohol-containing agents as they may cause protein fixation.

PREPARATION BEFORE CLEANING

If possible, disassemble the product into its individual components. Inspect for damage and wear.

CLEANING & DISINFECTION: BY MACHINERY (RECOMMENDED)

Only perfectly maintained, checked, calibrated and released cleaning and disinfection equipment according to EN ISO 15883 with valid check mark (CE mark) may be used. When machine cleaning and disinfecting several different products/individual components, ensure that the cleaning device is not overfilled and that the products/individual components do not contact each other. We recommend completing the cleaning and disinfection cycle according to the specifications of the manufacturer of the device and/or the cleaning agent. Make sure that the cleaning agent is compatible with the products and disinfectants that might be used. The pH of the cleaning agent and disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogen hydrocarbons. We recommend the following products: Cleaning agent: Dr. Weigert, neodisher MA/Neutralizing agent: Dr. Weigert, neodisher Z. The instructions for use of the manufacturer of the cleaning agent and disinfectant must strictly be observed. In particular, note the concentration and contact times.

Description of the recommended cleaning process:

Pre-cleaning: The pre-cleaning should be completed in a water bath (at least drinking water quality) without other cleaning additions using a cleaning brush below the water surface to achieve both sufficient cleaning of the products/individual components while avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product/individual component; however, it should be no less than 1 minute.

Cleaning tool for pre-cleaning: Flexbrush cleaning brush, REF 605 254/Mirasuc brush, REF 605 216

Information on the cleaning and disinfection device: Type: Miele G7881

Description of the recommended program:

Cleaning program:	VarioTD
Cleaning temperature:	50.2 - 56.4°C
Hold time:	10:54 min
Disinfection temperature:	65.1 - 94.3°C
Hold time:	09:57 min
Flushing technology:	Fresh water system/demineralized water
Cleaning agent:	Dr. Weigert, neodisher MA, Dosing: 30 mg (3g/l)
Neutralizer:	Dr. Weigert, neodisher Z, Dosing: 1 - 2 ml/l
Aids:	Insert rack for impression trays

The water should be of drinking water quality. Flushing and drying should be completed by machines. The machine cleaning should not exceed a temperature of 60°C. Ensure that visible contamination has been completely removed. If visible contamination remains after cleaning and disinfection, repeat these steps. Before further processing, the product must be free from any residues and dry.

CLEANING: MANUAL

The manual cleaning should be completed in a water bath (at least drinking water quality) with the cleaning agent specified below using a brush below the water surface to achieve both sufficient cleaning of the products/individual components while avoiding fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product/individual component; however, it should be no less than 1 minute. Afterwards, rinse the products/individual components with water (at least drinking water quality) for at least 20 seconds.

Description of the recommended cleaning process:

Pre-cleaning:	Remove surface contamination using a disposable cloth
Cleaning:	Insert the instrument into the strainer bowl. Hang the strainer bowl into the cleaning bath with cleaning solution.
Cleaning agent:	ID 213, company Dürr Dental, concentration: 2 %
Cleaning tool:	Flexbrush cleaning brush, REF 605 254/Mirasuc brush, REF 605 216
Rinsing:	Rinse with water (at least drinking water quality) for 20 sec.
Drying:	Dry at room temperature

The manual cleaning should not exceed a temperature of 45°C. Ensure that visible contamination has been com-

pletely removed. If visible contamination remains after the cleaning process, cleaning must be repeated.

DISINFECTION: MANUAL

For the manual disinfection, please use only approved disinfectants with tested effectiveness (CE mark, VAH-/DGHM-listed). Put the products/individual components into corresponding disinfection baths according to the disinfectant manufacturer's specifications. Ensure that the products are sufficiently covered and do not contact each other. The instructions for use of the disinfectant manufacturer must be strictly observed. In particular, observe the concentration and contact times. Application with the following agent is recommended: Dürr Dental ID 213 instrument disinfection, concentration 2%

Description of the recommended disinfection process:

Completion of proper cleaning preparation and manual cleaning

Disinfection: Insert the instrument into the strainer bowl. Hang the strainer bowl in the disinfection bath with disinfection solution.

Disinfectant: ID 213, company Dürr Dental, concentration: 2 %

Cleaning tool: Flexbrush cleaning brush, REF 605 254/Mirasuc brush, REF 605 216

Contact time: 5 min

Rinsing: Rinse with water (at least drinking water quality) for 20 sec.

Drying: Dry at room temperature

Make sure that the disinfectant is compatible with the products and cleaning agents that might be used. The pH of the disinfectant should be between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogen hydrocarbons. Please observe thermal restrictions that might apply, taking into account the disinfectant manufacturer's specifications. Afterwards, rinse the products/individual components with water (minimum drinking water quality) for at least 20 seconds. Before further processing, the product must be free from any residues and dry.

DRYING

If you process the products manually, the products/individual components may be dried by blowing them off with filtered, oil-free compressed air according to DIN ISO 8573-1 (medical cleanliness class) or at room temperature. The drying in machine processing should be directly completed by the devices using hot air. With machine drying, the temperature should not exceed 100°C.

MAINTENANCE, CONTROL AND CHECK

There are no special maintenance requirements for the product. Re-assemble disassembled products/individual components. Always perform a visual inspection for contamination, damage, wear and deformation before and after the individual work steps. Damaged or corroded products must no longer be used. If the product/individual component is not visibly clean, the entire processing must be repeated or the product/individual component must be disposed of properly.

PACKAGING

Pack and seal the product in suitable and standardized transparent sterile packaging (sterilization bags). Observe the instructions of the sterilization bag and sealing machine manufacturers as well as the current, normative requirements. Products/individual parts that have not been sterilized in bags must be used immediately.

STERILIZATION

Only tested steam vacuum autoclaves may be used. When sterilizing several different products/individual components, ensure that the autoclave is not overfilled and that the products/individual components do not contact each other. The following sterilization cycles may be completed: Steam sterilization, 134°C, hold time 5 minutes or steam sterilization, 121°C, hold time 20 minutes. To dry the products/individual components, the drying cycle of the autoclave should be set. Use the autoclave manufacturer's instructions for use.

STORAGE

To maintain the sterility, the products should be stored in standardized sterilization bags in a dry, clean place until they are used. If the sterile packaging is damaged, the products must be processed again before use.

ADDITIONAL INFORMATION

Mark the sterilized products according to the legal and national regulations. The recommended storage duration for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling. The user must ensure within the scope of their quality management system that specified processing cycles that might apply (see Limitations in the processing) are not exceeded. The manufacturer and the competent authority of your member state shall be immediately notified all serious incidents occurring in connection with the product. Recommendations and regulations for disposal of the products vary by region. Ask your competent disposal company for the current regulations in your region.

CONTACT TO THE MANUFACTURER

Phone: (+49) 203/99269-0, Fax: (+49) 203/299283, Email: info@hagerwerken.de

The instructions specified above are validated by the medical device manufacturer as suitable for the preparation of a medical device for its reuse. The processor is responsible for ensuring that the actually completed processing with the equipment used, materials and personnel in the preparation facility achieves the desired result. To this end, verification and/or validation and routine monitoring of the procedure are required.

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