Master all your challenges

initiaL REVO MC

EN Instructions For Use

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APRIL 2026

Description

Dental Ceramic

Intended purpose

GC Initial REVO MC is a veneering ceramic with EN ISO 6872 for the ceramic veneering of precious and non-precious metal substructures with a coefficient of thermal expansion of the metal substructure in the range of $13,8 - 14,9 \cdot 10^{-6} \text{ K}^{-1}$ (25-500°C). Also suitable for the sintering technique to create inlays, onlays and veneers on refractory dies.

Indications for use

For characterisation and veneering of dental alloy frameworks:

- Crown and Bridgework (incl reproduction of gum tissue)
- Implant superstructures

For the production of reconstructions without a framework:

- Jacket Crowns, Inlays, Onlays and Veneers

Clinical Benefits

- Restoring and/or improving masticatory function
- Improving teeth aesthetics

Contraindications

- For patients with parafunctions (e.g. bruxism)
- Multi-unit framework free restorations
- Veneering of alloys in the non-approved CTE range (CTE= coefficient of thermal expansion of the ceramic veneering material or the used alloys. The respective CTE must be coordinated with one another)
- It may be contra-indicated in patients with known allergies to specific metal alloys

Intended users / Special training

- Dental technicians
- Dentists (clinical procedure)
- No special training required

Patient target group

Patients with permanent teeth

Composition

GC Initial REVO MC

Powders:

SiO₂ 55-65%, Al₂O₃ 12-16%, K₂O 12-16%, Na₂O 3-5%, Li₂O <1%, MgO <1%, CaO 1-3%, BaO 1-2%, B₂O₃ 1-2%, TiO₂ <1%, ZrO₂ <1%, P₂O₅ <1%, CeO₂-CeF₃ <1%, SnO₂ <1%, Pigments 0.1-3.0%

Modelling Liquid:

H₂O + Zinc chloride, Butanediol

GC Initial Spectrum Glaze Powder

SiO₂ 61-65%, Al₂O₃ 8-10%, K₂O 5-7%, Na₂O 7-10%, CaO 1-2%, BaO 2-3%, B₂O₃ 7-10%, TiO₂ <1%, ZrO₂ <1%, P₂O₅ 1-2%, CeO₂-CeF₃ <1%, Pigments 0.1-1.0%

GC Initial Spectrum Glaze LiquidH₂O + Propylene glycol**GC Initial MC Shoulder Liquid**H₂O, Zinc chloride, Tylose**Performance characteristics**

The product's performance characteristics satisfy the requirements of the intended use.

Direction for use**1. Preparing the substructure**


- Use non-precious, semi-precious, or precious alloys with a coefficient of thermal expansion: $13.8 - 14.9 \times 10^{-6} \text{ } ^\circ\text{C}^{-1}$.
- Shape the framework to mirror the tooth's contour in a reduced form.
This approach ensures consistent layer thickness and uniform support for the veneering ceramic.
- Follow the alloy manufacturer's instructions for fabricating and oxidizing the framework.

2. Masking of the substructure

- **Bonder application (optional)**
 - GC INmetalbond is used as a thin layer between the alloy and the first Opaque layer. It blocks escaping metal oxides and neutralizes differences in the expansion coefficient (CTE). INmetalbond can be used on all alloys within the given CTE range.
 - Apply the bonder in a thin layer but masking the framework completely. The fired bonder should have a yellowish, slightly shiny surface.
 - For more information check the dedicated Instructions for Use
- **Opaque application**
 - **First Opaque Layer:**
 - Apply a thin, even layer completely masking the coping.
 - After firing, the first opaque will appear shiny with partial coverage of the applied surface.
 - **Second Opaque Layer:**
 - Apply a thin, even layer over the first opaque, ensuring complete masking.
 - After firing, the opaque will have a slight shine.

3. Body & Enamel Build Up

Build the Body to full contour and reduce for incisal layer.

**4. Shoulder Application**

- Make sure to isolate the die in the shoulder area.
- Reduce the metal frame in the shoulder area to a minimum. The shoulder is covered up to the preparation border with the Shoulder Universal material mixed with the corresponding Body mixture selected for the tooth shade.

5. Glaze Firing

- For glazing use the GC Initial Spectrum Glaze and Glaze Liquid.
- For enhanced individualization, GC Initial IQ Lustre Pastes ONE/NF Gum and GC Initial Spectrum Stains can be applied to achieve highly aesthetic, lifelike results – simply by painting.

Note:

For more tips and tricks please consult the GC Initial REVO MC Technical Manual on our dedicated website www.gc.dental/europe.

Firing advice

FIRING CHART REVO MC						
	Preh. Temp.	Drying Time	Raise of Temp.	V	Final Temp.	Holding Time
INmetalbond	550°C	6 min.	80°C/min.	+	980°C	1 min.
REVO_/P/O_1. Firing (wash)	550°C	6 min.	75°C/min.	+	980°C	1 min.
REVO_/P/O_2. Firing	550°C	6 min.	75°C/min.	+	950°C	1 min.
REVO_Shoulder	550°C	2 min.	75°C/min.	+	945°C	1 min.
REVO_BODY_1. Firing	550°C	6 min.	55°C/min.	+	915°C	1 min.
REVO_BODY_2. Firing	550°C	6 min.	55°C/min.	+	905°C	1 min.
Glaze Firing (without Stain/Glaze)	600°C	2 min.	55°C/min.	–	890-920°C*	1 min.
Glaze Firing (with Stain/Glaze)	550°C	2 min.	55°C/min.	+	870-910°C*	1 min.
Individualisation with Lustre Paste ONE/Lustre Paste NF Gum	480°C	2 min.	45°C/min.	+	810°C	1 min.

Note:

The given firing temperatures are recommendations only and may differ depending on the furnace. The temperatures may need to be adjusted to the furnace you are using and according to the result after firing. Always make sure your furnace is correctly calibrated.

* The firing temperature is determined by the furnace used and the intended gloss level.

Disposal

Dispose of all waste according to local regulations.

Storage

Recommended for optimal performance, store at 4-25°C (39-77°F).

Side effects

There are no known side effects to date.

Residual risks

- Technical complications: chipping
- Biological complications: secondary caries

Precaution, warnings and limitation of use

Inhalation of ceramic dust released during processing must be avoided by using suitable suction devices and, if necessary, protective masks. Personal Protective Equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the Safety Data Sheets available at: <http://www.gc.dental/europe> or for the Americas: <http://www.gc.dental/america>. They can also be obtained from your supplier.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at klema.qm@klema.at (Basic UDI-DI: (01)9010782MD001CE)

Package**GC Initial REVO MC STARTER SET**

3 Paste Opaque: Light / Medium / Dark
2 Body: Body A Light / Body A Dark
3 Enamel: Light / Medium / Dark
3 Effect: Opal Booster / CL-F / TN
1 Shoulder Universal
1 Modelling Liquid
1 MC Shoulder Liquid
1 Spectrum Glaze Powder GL
1 Spectrum Glaze Liquid

GC Initial REVO MC REFILLS

Paste Opaque: Light / Medium / Dark / Gum Universal / Bleach
Body: Body A Light & Dark / Body B Light & Dark / Body C Light & Dark / Body D Light & Dark
Enamel: Light / Medium / Dark
Effect: TN / Opal Booster / CL-F
Shoulder Universal
Gum: Light / Dark
Bleach: Body / Enamel
Modelling Liquid

Reporting Obligation according to EU MDR

In the case of serious incidents related to the device, please report it to KLEMA Dentalprodukte GmbH. (Klema.qm@klema.at) and to the competent authority of the Member State in which the user/patient resides.

Note:

The current Instructions for Use are available in the download section of the <https://www.gc.dental/> website. To access the appropriate version, please select your region from the footer.

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EU: Klema Dentalprodukte GmbH
ALSIP, IL 60803 U.S.A.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
GC ASIA DENTAL PTE. LTD.
GC AUSTRALASIA DENTAL PTY.LTD