KAVITAN™ LC

Light-Cured Glass-Ionomer Cement

en

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Contraindicaciones:
The product is not recommended for mechanically stressed Class IV restorations. Class II restorations with conventional cavity preparation in permanent teeth, and for restoration of occlusal cusps. Do not use the product in case of known hypersensitivity to any component of the product.

INSTRUCTION FOR USE:

Powder : liquid mixing ratio	2.8 : 1 (g:g)
Mixing time	30 - 45 sec
Working time	min. 3 min
Photopolymerization - light curing time	20 sec

The above mentioned values apply to material processed at 23°C and air humidity of 50%.

Cement preparation:

- Measuring 2.8 : 1 (g : g), i.e. 1 measuring scoop of powder to 2 drops of liquid. Optimal mixing time is 30-45 seconds.
- Shake powder in the bottle before use. Do not compress powder in the measuring scoop. Overfill the measuring scoop with powder and then remove surplus powder using the edge of the cutting-off insert in the bottle neck.
- Hold the bottle vertically when doing the liquid. You will get the optimum drop size when you squeeze the bottle gently. Any dried material must be removed from the dropper.
- Close the bottles immediately after use. The material should be protected from moisture.
- Pour the powder and liquid on a disposable smooth paper pad or a glass plate. Powder can be split in 2-4 parts. Mix powder with liquid using a plastic or stainless spatula. Mixing time should be as short as possible (30-45 sec at most). Properly prepared cement is homogeneous and plastic.

Principles of clinical application:
- Conventional preparation of carious lesion including macro-removal (do not prepare carious cavities and abrasions lesions, only clean the cavity surface).

- For pulp capping use calcium hydroxide base material. Do not use eugenol-based cements, that inhibit polymerization.

- Improved adhesion can be achieved by application of conditioner - in accordance with manufacturer's instructions for use.

- Pour the mixed cement directly into the instrument with high-quality surface, or with special applicators (Centrix tips, etc.).

- Setting of the material under a transparent matrix can improve a surface quality of the restoration.

- Light cure for 20 sec. (470 nm wavelength) in layers with a maximum thickness of 2 mm.

- Excess of set material can be trimmed carefully immediately after polymerization using sharp rotary instruments or a suitable hand instruments.

- Sealing of the restoration with protective varnish (LC VARNISH) improves the surface quality.

Caution:
The material adheres to instruments and it must be removed before setting.

Wash:
Liquid: H319 Causes severe eye irritation. H315 Causes skin irritation. H317 May cause an allergic skin reaction.

P262 Do not get in eyes, on skin, or on clothing. P270 Do not eat, drink or smoke when using this product.

Waste material disposal:
Residues of the unused material as well as the packaging material can be disposed of together with ordinary waste material.

Storage:
Store in a dry and dark place at 5 °C to 25 °C, in a well sealed container.

Shelf life:
2 years

Indication on the package: Exp YYYY-MM (YYYY = year, MM = month).

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