SpofaDental

A Kerr Company

 **SpofaDental a.s.**, Markova 238, 506 01 Jičín, Czech Republic
+ 420 493 583 111, www.spodadental.com. Made in Czech Republic

Silisan[®] N

Special liquid for phosphate-bonded investment materials

Characteristics:

Silisan N is a special liquid improving the properties of phosphate-bonded investment materials.

Application:

The investment material with high strength and the optimum setting and thermal expansion will be created by mixing with the powdery investment material Silikan. The concentration of Silisan N is adjusted according to the cast alloy type. The recommended water solutions of Silisan N are specified in the instruction leaflets of the individual investment materials. Silisan N shake well before use.

Warning:

When using metal casting moulds (rings), it is necessary to line the inside of the mould with the ceramic strip (KP-50) which allows the investment material expansion, eliminating thus the cast alloy contraction. EUH210 - Safety data sheet available on request.

Waste disposal:

The residue of unused product and contaminated packaging material should be transferred to a company authorized to dispose of waste materials. Observe your local legal regulations.

Storage:

Store in a dry place at a temperature of 20 ± 5 °C, in a well-closed package. There is a danger of preparation degradation at lower temperatures (below 5 °C).

Expiry period: 2 years

Package: 1000 g

Emergency telephone number:

United Kingdom, National Poisons Information Service (Newcastle Unit)

Claremont Place Newcastle-upon-Tyne, Newcastle, +44 191 2606182/+44 191 2606180 24H

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