

initial™ REVO MC

The [R]EVOLutionary Metal-Ceramic System for Modern Dentistry

Maximise your efficiency. Optimise your outcomes.



Since 1921
100 years of Quality in Dental



GC Initial™ REVO MC


The [R]EVOLutionary Metal-Ceramic System

Dentistry is evolving – and so are your needs.

With the rise of monolithic and minimally veneered materials, workflows must be faster, more efficient and predictable without sacrificing aesthetics. Initial REVO MC defines a new standard in metal-ceramic dentistry. It combines the proven strength of a metal framework with a contemporary, minimal ceramic design.

The result? Reduced veneering layers, streamlined processes, and exceptional aesthetics.

Benefits at a glance



Efficiency, redefined

Simplified workflows with fewer steps
Greater productivity and fewer remakes



Beauty without complexity

Stunning results with minimal layers
All shades, including bleach and gum tones



Reliable performance

Superior strength and durability to stand the test of time
Highly predictable shade reproduction

This is not your usual metal-ceramic system.
A new benchmark - inspired by the efficiency, strength, robustness and reproducibility of modern all-ceramic solutions.




[R]EVOLutionary colour concept

Using only 3 Paste Opaques, 8 Body Powders, and 3 Enamels, the full A-D spectrum can be accurately matched through a controlled, layer-based approach.



Each time a spot-on colour

The Initial REVO MC concept delivers a streamlined, high precision system for reproducing all 16 A-D shades with minimal material complexity.



Standardised mixing recipes

The standardised mixing ratios of the Body Powders ensure consistent, reproducible results, easily achieved with the help of a ceramic measuring spoon, making the workflow both efficient and reliable.









Thickness-independent colour precision

These standardised samples with ceramic layers from 0.5 to 2.0 mm thick, show the reliable colour result, regardless of the thickness. An efficient process with a uniform outcome.

- High-end shade accuracy with fewer products
- Full A-D shade coverage plus bleach, gum shades, and effect powders




OPAQUE							
BODY							
ENAMEL							
EFFECT							
SHOULDER							
BLEACH							
GUM							


initial

[R]EVOLutionary composition

Engineered with an *enhanced crystal structure* for superior durability, this advanced composition combines increased strength with exceptional aesthetic appeal.




Initial REVO MC - dense leucite reinforced network
Conventional porcelain fused to metal (PFM) - glass ceramic




- High robustness fitting the non-precious alloy
- Natural light dynamics enhancing the aesthetic value

[R]EVOLutionary Opaquers



- Perfect coverage with reliable bonding
- Strong adhesion with thinner and fewer layers

[R]EVOLutionary Modeling Liquid




- Smooth, intuitive handling
- Minimal shrinkage and dense, homogeneous firing results

Endless Possibilities

Compatible with your favourite modifiers

Stay in your flow – this is an upgrade with absolutely no limits! Initial REVO MC works seamlessly with the immensely popular Initial Lustre Pastes ONE and Initial Spectrum Stains for ultimate individualisation possibilities. No compromises, no learning curve – just enhanced performance with the modifiers that make your work shine.



Initial Lustre Pastes ONE are ready-to-use pastes with a fine thixotropic property, for easy and precise application. They can be used in combination with Initial REVO MC, below them or on top.



Initial Spectrum Stains can also be mixed directly into the Initial REVO MC veneering ceramics for internal modification or staining. This combination opens up a broader artistic palette, elevating creative freedom and enabling truly personalised, high-impact results.

You can bring Initial REVO MC restorations even more to life with simple, precise painting effects using Initial Lustre Pastes ONE, Initial Lustre Pastes NF GUM and Initial Spectrum Stains. Achieve the highest level of individualisation in every brushstroke – with no complex layering or extra firing cycles required. Fast, easy and secure: stunning effects in a single firing!



The New Standard in Metal-Ceramic Dentistry

Initial REVO MC is more than a product – it's a concept that redefines simplicity without compromising beauty.

Maximise your efficiency. Optimise your outcomes.

**Ready to Experience the [R]EVOLution?
Contact us today to learn more or request a demo.**





Sustainable Packaging



Eco-friendly for a clean world... and clean conscience

The packing of Initial REVO MC is designed for a reduced environmental impact

- Fully recyclable materials
- Reusable jars
- Minimised packaging volume and weight

initial™ REVO MC

The [R]EVOLutionary Metal-Ceramic System for Modern Dentistry

GC Initial REVO MC STARTER SET

3 Paste Opaque: Light / Medium / Dark
2 Body: Body A Light / Body A Dark
3 Enamel: Light / Medium / Dark
3 Effect: Opal Booster / CL-F / TN
1 Shoulder Universal
1 Modelling Liquid
1 MC Shoulder Liquid
1 Spectrum Glaze Powder GL
1 Spectrum Glaze Liquid

GC Initial REVO MC REFILLS

Paste Opaque: Light / Medium / Dark / Gum Universal / Bleach
Body: Body A Light & Dark / Body B Light & Dark / Body C Light & Dark / Body D Light & Dark
Enamel: Light / Medium / Dark
Effect: TN / Opal Booster / CL-F
Shoulder Universal
Gum: Light / Dark
Bleach: Body / Enamel
Modelling Liquid



We gratefully acknowledge:

- ZTM Stefan Roozen - Austria
- MDT Mark Bladen - UK
- MDT Patric Freudenthal - Sweden

for sharing the images of their craftsmanship showed in this leaflet.

GC EUROPE N.V.

Head Office
Researchpark,
Haasrode-Leuven 1240
Interleuvenlaan 33, B-3001 Leuven
Tel. +32 16 74 10 00
Fax. +32 16 40 48 32
info.gce@gc.dental
<https://www.gc.dental/europe>

GC UNITED KINGDOM Ltd.

Coopers Court Newport Pagnell
Buckinghamshire MK16 8JS
United Kingdom
Tel. +44 1908 218 999
Fax. +44 1908 218 900
info.uk@gc.dental
<https://www.gc.dental/europe/en-GB>