

NO.4. 40

CT

ProphyCare[®]

Directa's products are developed and evaluated by our selected team of highly-qualified and renowned Swedish and international dentists. It is our goal to make life easier for the dentist in the daily clinical work.







The perfect starter, effective removal of stubborn stains

For discoloured surfaces and thick layers of plaque

Prod No Product Description

 731113
 ProphyCare Prophy Paste, blue RDA 250, 125 micron, tubes, 96 g 60 ml

 690123
 ProphyCare Prophy Paste Single Dose, blue RDA 250, 144x 2 g

 690119
 ProphyCare Prophy Paste Single Dose, mix, 144x 2 g



Less abrasive, but still effective

For thick layers of plaque and discoloured surfaces

Prod No Product Description

731114 ProphyCare Prophy Paste, green RDA 170	, 70 micron, tubes, 97 g 60 ml
690122 ProphyCare Prophy Paste Single Dose, gre	en RDA 170, 144x 2 g
690119 ProphyCare Prophy Paste Single Dose, mix	x, 144x 2 g





Gentle cleaning, tissue friendly stain remover

For small amount of plaque, root surfaces and lightly discoloured teeth

Prod No Product Description

731112 ProphyCare Prophy Paste, red RDA 120, 20 micron, tubes, 84 g 60 ml 690121 ProphyCare Prophy Paste Single Dose, red RDA 120, 144x 2 g

690119 ProphyCare Prophy Paste Single Dose, mix, 144x 2 g

Profylaxpasta CCS Finuts Finuts Frinuts For final polishing For final polishing Four le polissage final Pour le polissage final Pour le polissage final Pour le polissage final



Ideal for final polishing

For final polishing of teeth and fillings

Prod No	Product Description
731111	ProphyCare Prophy Paste, yellow RDA 40, 2 micron, tubes, 84 g 60 ml
690120	ProphyCare Prophy Paste Single Dose, yellow RDA 40, 144x 2 g
690119	ProphyCare Prophy Paste Single Dose, mix, 144x 2 g

Maximum hygiene

Always a fresh sealed package No risk of cross contamination, and the highest hygiene guarantee for your patient.

Adhesive base

The single dose cup has an adhesive base that easily sticks to the glove or any other surface.

No waste

The shallow packaging with rounded edges is designed to minimize waste. All the content can easily be collected by the polishing cup.

Perfect size

Each package contains the ideal quantity for one treatment.

ProphyCare[®] Prophy Paste Single Dose

Prophy Paste single dose is available in assortment or refill packages. The selection guide included in each box will help you to choose the appropriate value for the treatment. The box is designed to keep the single dose cups tidily in place and quickly give an overview of its content allowing users to select what is needed for next treatment. Choose Prophy Paste after evaluating for degree of discoloration, exposed dentin or cementum, aesthetic and other restorative materials, areas of demineralisation etc. Use at least two different Prophy Paste RDA values for each treatment - RDA 250 and/or RDA 170 for cleaning and RDA 120 and/or RDA 40 for polishing.

Prod No Product Description

690123	ProphyCare Prophy Paste Single Dose, blue RDA 250, 144x 2 g
690122	ProphyCare Prophy Paste Single Dose, green RDA 170, 144x 2 g
690121	ProphyCare Prophy Paste Single Dose, red RDA 120, 144x 2 g
690120	ProphyCare Prophy Paste Single Dose, yellow RDA 40, 144x 2 g
690119	ProphyCare Prophy Paste Single Dose, mix, 144x 2 g
690125	ProphyCare Prophy Paste PRO Single Dose, 144x 2 g





Selection Guide

	Blue RDA 250	Green RDA 170	Red RDA 120	Yellow RDA 40
Starter				
Very discoloured surfaces				
Thick plaque				
Discoloured teeth				
Small amounts of plaque				
Root surfaces				
Slightly discoloured teeth				
Fillings				
Pediatric				
Cervical toothbrush abrasion lesions				•
Implants				•
Final polishing				
After air polishing treatment				
Final polish				



ProphyCare[®] Prophy Paste Blue Coarse, Paraben Free RDA 250, 125 micron

Perfect starter For discoloured surfaces and thick layers of plaque Effective removal of stubborn stains caused by e.g. coffee, tobacco Abrasive agent: pumice Contains 0,1% fluoride Use in combination with other ProphyCare Prophy Paste

 Prod No
 Product Description

 731116
 ProphyCare Prophy Paste, blue, Paraben Free RDA 250, tubes, 96 g 60 ml



ächen und besonders starken



RDA 250



ProphyCare[®] Prophy Paste Red Fine, Paraben Free RDA 120, 20 micron

For removing plaque and surface discolourations For small amount of plaque, root surfaces and lightly discoloured teeth Excellent for polishing Abrasive agent: silica Contains 0,1% fluoride Use in combination with other ProphyCare Prophy Paste

 Prod No
 Product Description

 731115
 ProphyCare Prophy Paste, red, Paraben Free RDA 120, tubes, 84 g 60 ml



1.11111112019-05111111

Low abrasion and still excellent cleaning effect

Prophy Paste PRO low RDA value of 20 offers effective stain removal properties at the beginning of the treatment and decreases progressively to a lower RDA value of 5, giving the paste its final polishing properties at the end of the treatment.

One step procedure

The same paste is used for stain removal and final polishing. ProphyPaste PRO is the one-step prophy paste for outstanding cleaning & polishing results that offers decreased treatment time.

Remineralizing agent

Contains fluoride to remineralize enamel after polishing, thus supporting caries prevention.

Integrated abrasion ariability

Prophy Paste PRO's high cleaning power converts into gentle polishing as the particles fragment.



PIRECTA

ProphyCare[®] Prophy Paste PRO, Lilac Universal

Prophy Paste PRO is the one-step prophy paste for outstanding cleaning and polishing results that offers decreased treatment time. The same paste is used for stain removal and final polishing. Prophy Paste PRO's high cleaning power converts into gentle polishing as the particles are fragmented. Prophy Paste PRO's RDA value of 20 offers effective stain removal properties at the beginning of the treatment and decreases progressively to a lower RDA value of 5, giving the paste its final polishing properties at the end of the treatment. Prophy Paste's mild mint taste and appealing texture is pleasant for patients of all ages and does not splatter, liquify or separate and is easy to rinse-off. Contains fluoride to remineralize enamel after polishing, thus supporting caries prevention. Available in single doses and tubes.





ProphyCare[®] Prophy Powder Effective powder for stain and plaque removal

ProphyCare Prophy Powder is the new effective powder for removal of stains and plaque. The Prophy Powder effectively cleans the tooth surface, penetrating into deep pits, fissures and interproximal areas. Heavy stains are thoroughly and gently removed. ProphyCare Prophy Powder contains anti-clogging agents for a trouble-free performance in any air polishing unit, ensuring an even and maximum flow of powder is continuously applied on each patient. Safe for use in all air-polishing units.

Using Prophy Powder is time saving, safe and effective for stain and plaque removal thus minimizing the operator and patient fatigue. Its lemon aroma eliminates any salty aftertaste leaving the mouth feeling refreshed.

 Prod No
 Product Description

 741000
 ProphyCare Prophy Powder, lemon taste, supragingival, 1 x 250 g

ProphyCare[®]

NyCare[®]

^{yPowder}

Ingival Dicarbonate

190504

2021-09

Vasby, Sweden

CETi

Accessible & effective

Safe and effective way to remove stains. ProphyCare ProphyPowder will clean the tooth surface, penetrating into deep pits, fissures and interproximal areas and make you save time.

Removes heavy stains

Removes stubborn and heavy stains, with minimal damage to teeth.

Fresh lemon taste

Its lemon aroma eliminates any salty aftertaste leaving a fresh feeling in the mouth.

Anti-clogging

Contains anti-clogging agents for maximum flow and troublefree performance in your air polishing units.

Interior webbing

Inner structure design minimizes splatter and increases friction for efficient stain removal.

Exterior ridges

Outer ridges optimize interproximal cleaning and splatter control. Soft flare improves subgingival access.

Soft edge

Flexible flare for improved subgingival access and efficient plaque removal.

Protective base

Wider base designed to prevent prophy paste from reaching the contra angle.



ProphyCare[®] Prophy Cup Optimal cleaning, superior paste holding

The exclusive soft and flexible ProphyCare Prophy Cup flares and flexes to contours of teeth to aid in selective stain removal and interproximal cleaning with less pressure and operator fatigue. Optimized size for easy access in hard to reach areas without compromising the paste holding capacity. The unique design is suitable for treatment on adults and pediatric use and to substitute pointed polishers.

Prod No	Product Description
690131	ProphyCare Prophy Cup

690131ProphyCare Prophy Cups, screw-in, 144 pcs690132ProphyCare Prophy Cups, latch, 144 pcs



ProphyCare[®] Prophy Angle Efficient cleaning, enhanced comfort

Disposable prophy angle in polycarbonate with soft, flexible 100% latex-free cup for low speed prophy attachment handpieces. Prophy Angle minimize prophy paste particles to enter the handpiece and damage the mechanism. Smooth running, silent, vibration-free and no overheating. With its ergonomic shape and slim neck design ProphyCare Prophy Angle enhances access to posterior teeth and increases comfort for the patient while reducing hand fatigue and repetitive wrist motion for the hygienist.

Prod NoProduct Description690130ProphyCare Prophy Angle, 144 pcs

Interior webbing

Inner structure design minimizes splatter and increases friction for efficient stain removal.

Exterior ridges

Outer ridges optimize interproximal cleaning and splatter control. Soft flare improves subgingival access.

Maximum hygiene

Designed to minimize risk of cross contamination and reduce maintenance of the handpiece.

Ergonomic design

Slim neck, ergonomic angle and round head for enhanced visibility and comfortable access to all areas of the mouth.



ProphyCare[®]

Design By Dentists™

Directa's products are developed and evaluated by our selected team of highlyqualified and renowned Swedish and international dental professionals with the aim of making life easier for the dental professionals in the daily clinical work.

DIRECTA directadental.com

DIRECTA AB P.O. Box 723, 194 27 Upplands Väsby, SWEDEN +46 8 506 505 75

DIRECTA INC 64 Barnabas Road, Unit 3 Newtown, CT 06470, United States 203.491.2273

2020-01-0

0

DIRECTA GERMANY Hammergutstr. 11 92245 Kümmersbruck, GERMANY

info@directadental.com www.directadental.com www.directadentalgroup.com www.facebook.com/directadental www.instagram.com/directadental www.youtube.com/user/DirectaAB

