

KENDA C.G.I.

3-STEP POLISHING SYSTEM

- Comprehensive polishing system with three grits
- Covering the whole process from contouring to high shine finish
- Flexible, adapting to tooth shape
- White polisher produces extra cutting efficacy for quick removal of excess material



SET REF 900.ASS.012

16



REF	reorder	905-C	908-C	909-C	910-C	905-M	908-M	909-M	910-M	905-F	908-F	909-F	910-F
ISO	1: 638 204 2: 638 204... 3: 638 204...	030 533 060	243 533 050	303 533 100	030 533 090	030 523 060	243 523 050	303 523 100	030 523 090	030 493 060	243 493 050	303 493 100	030 493 090
	∅ × L mm	6,5 x 10	5 x 10	10 x 1	9 x 8,5	6,5 x 10	5 x 10	10 x 1	9 x 8,5	6,5 x 10	5 x 10	10 x 1	9 x 8,5
	μ grit												
	color												
	rec.												
	max.												
	μ grit												
	color												
	rec.												
	max.												
	Refill Box of:	N Pressure:				Autoclavable:				Shank:		Medical device:	
	6, 12, 25, 100	0,3 - 2N				135°C, 275°F				RA, stainless steel		Class IIa	

